



- Sinolinks is a high-tech medical manufacturer in mechanical suturing instruments, surgical repairing materials, polymer ligating clips, and minimally invasive surgery related products in the field of gastro-enterology, colo-rectal surgery, general thoracic surgery, urology and gynecology.
- Company has top-ranking talent team and modern production & testing facilities, we concentrate various advantage resources to ensure high quality products and services.

Vision

- Exert ourselves to improve medical technology, and delicate ourselves to surgery in order to offer comprehensive range of solutions.
- Long for products with the highest quality persistently, pursuit the harmony between enterprise and clients, and spear no effort to exceed customers' requirements and expectations.
- Have the courage to undertake social responsibility, and create value for employees, customers, and shareholders.

Honor

- Forbes China – 2015 China top 100 unlisted potential enterprises
- China Finance Summit – 2015 the most potential developing enterprises
- Asia Brand Ceremony – Asian famous brands awards enterprise, General Manager was awarded “top 10 China (Industry) brand innovation characters”
- National high new tech enterprises
- Member of the national standardization technical committee (SAC/TC94) surgical instrument
- China FDA, CE 1097, ISO13485, FDA (510k)
- Products are widely used in more than 1000 domestic and foreign hospitals

Sinolinks Medical Innovation, Inc.

Add: NO.10 Ziwei Road, Zhonglou Zone, Changzhou, Jiangsu, China. 213032
Tel: +86-519-88021691
Fax: +86-519-88021698
Website: www.sinolinks.com.cn
E-mail: sales@sinolinks.com.cn



Science and Technology for Human Health

SINOLINKS MEDICAL INNOVATION, INC.



Disposable Circular Stapler
DCS A/E/H

Code	O.D (mm)	I.D (mm)	Staple	Height (mm)
DCS H 17	17	10	16	4.5
DCS H 19	19	12	18	4.5
DCS A/E/H 21	21	13	18	4.8
DCS A/E/H 23	23	15	18	4.8
DCS A/E/H 25	25	17	20	4.8/5.2
DCS A/E/H 27	27	18.5	22	4.8
DCS A/E/H 29	29	20	24	5.2
DCS A/E/H 31	31	22	28	5.2
DCS A/E/H 33	33	24	30	5.2

Disposable Hemorrhoidal Stapler

DHS B



Code	O.D(mm)	I.D(mm)	Staple	Height (mm)
DHS B-32	32	24	32	3.8
DHS B-34	34	26	34	3.8

Disposable Linear Stapler

Code	Colour	Length (mm)	Staple	Height (mm)
DLS B-30	Blue	30	11	3.8
	Green	30	11	4.8
DLS B-45	Blue	45	15	3.8
	Green	45	15	4.8
DLS B-60	Blue	60	21	3.8
	Green	60	21	4.8
DLS B-90	Blue	90	33	3.8
	Green	90	33	4.8

Linear Reload

Code	Colour	Height (mm)
DLR 30-B	Blue	3.8
DLR 30-G	Green	4.8
DLR 45-B	Blue	3.8
DLR 45-G	Green	4.8
DLR 60-B	Blue	3.8
DLR 60-G	Green	4.8
DLR 90-B	Blue	3.8
DLR 90-G	Green	4.8

Disposable Linear Stapler
DLS B



Code	Length (mm)	Staple	Height (mm)
DCC B-40	40	40	4.8

Disposable Curved Cutter Stapler
DCC B

Disposable Linear Cutter Stapler
DLC B



Disposable Linear Cutter Stapler

Code	Colour	Length (mm)	Staple	Height (mm)
DLC B-60	Blue	61	60	3.8
	Green	61	60	4.8
DLC B-80	Blue	81	84	3.8
	Green	81	84	4.8
DLC B-100	Blue	101	104	3.8
	Green	101	104	4.8

Linear Cutter Stapler Reload

Code	Colour	Height (mm)
DCR 60-B	Blue	3.8
DCR 60-G	Green	4.8
DCR 80-B	Blue	3.8
DCR 80-G	Green	4.8
DCR 100-B	Blue	3.8
DCR 100-G	Green	4.8

Disposable Endoscopic Linear Cutter

DEC



Disposable Endoscopic Linear Cutter

Code	Shaft Length (mm)	Suitable Cartridges
DEC A/B 60	60	DER A/B 30, 45, 60
DEC A/B 160	160	DER A/B 30, 45, 60
DEC A/B 260	260	DER A/B 30, 45, 60

Endoscopic Linear Cutter Cartridges DER


Code	Color	Suture Length (mm)	Staples (No.)	Height (mm)
DER A 30 V	Gray	31	48	2.0/2.0/2.0
DER A 30 L	White	31	48	2.5/2.5/2.5
DER A 30 H	Blue	31	48	3.5/3.5/3.5
DER A 45 V	Gray	46	66	2.0/2.0/2.0
DER A 45 L	White	46	66	2.5/2.5/2.5
DER A 45 H	Blue	46	66	3.5/3.5/3.5
DER A 45 M	Green	46	66	4.8/4.8/4.8
DER A 60 L	White	61	90	2.5/2.5/2.5
DER A 60 H	Blue	61	90	3.5/3.5/3.5
DER A 60 M	Green	61	90	4.8/4.8/4.8
DER B 30 G	Gray	31	48	2.0/2.0/2.0
DER B 30 T	Tan	31	48	2.0/2.5/3.0
DER B 30 P	Purple	31	48	3.0/3.5/4.0
DER B 45 G	Gray	46	66	2.0/2.0/2.0
DER B 45 T	Tan	46	66	2.0/2.5/3.0
DER B 45 P	Purple	46	66	3.0/3.5/4.0
DER B 45 B	Black	46	66	4.0/4.5/5.0
DER B 60 T	Tan	61	90	2.0/2.5/3.0
DER B 60 P	Purple	61	90	3.0/3.5/4.0
DER B 60 B	Black	61	90	4.0/4.5/5.0

Code	Working Length (mm)
VN1.5-120	120
VN1.5-150	150

Veress Needle
VN




Disposable Trocar
BVT



Code	Diameter (mm)	Working length (mm)
A TYPE Safety Tip Trocar	BVT A5	5
	BVT A10	10
	BVT A12	12
	BVT A15	15
B TYPE Pyramid Tip Trocar	BVT B5	5
	BVT B10	10
	BVT B12	12
BVT B15	15	100,150

Polymer Ligating Clips



Code	Size	Length of the clip (mm)	Scope of the closure Diameter (mm)
LC-1	ML	10	3-10
LC-2	L	13	5-13
LC-3	XL	17	7-16

Laparoscopic Applier

LC	Code	Size	Color
	CA B 5-S	Φ5mm x 330mm	Green
	CA B 10-M	Φ10mm x 330mm	Purple
	CA B 10-L	Φ10mm x 330mm	Golden

Disposable Fascial Closure System

DFCS

Code	Diameter (mm)	Total Length (mm)
DFCS -13	13	165



Code	Diameter (mm)	Volume (ml)
DSRB A-40	40	40
DSRB A-60	60	75
DSRB B-80	80	200
DSRB B-100	100	350
DSRB B-130	130	700



Disposable Specimen Retrieval Bag

DSRB

Disposable Hand Assist Tunnels (Incision Protectors)

Code	O.D (mm)	I.D (mm)	Channel Diameter (mm)	Length (mm)
40/30-15/25	40	30	15	25
50/40-15/20	50	40	15	20
50/40-25/25	50	40	25	25
50/50-25/30	50	50	25	30
SQB 60/50-25/25	60	50	25	25
(Silicon) 60/50-25/50	60	50	25	50
70/60-35/50	70	60	35	50
70/70-35/25	70	70	35	25
70/70-35/40	70	70	35	40

Code	O.D (mm)	I.D (mm)	Channel Diameter (mm)	Length (mm)
60/70-60/150	60	70	60	150
80/90-80/150	80	90	80	150
120/130-120/250	120	130	120	250
SQA 150/160-150/250	150	160	150	250
(PU) 180/190-180/250	180	190	180	250
220/230-220/250	220	230	220	250
270/280-270/250	270	280	270	250
320/330-320/250	320	330	320	250



Disposable Skin Stapler

DSS

Code	Staple Size (mm)	Staples (No.)
DSS A/B-15W	7.0x4.5	15
DSS A/B-15R	5.5x3.5	15
DSS A/B-25W	7.0x4.5	25
DSS A/B-25R	5.5x3.5	25
DSS A/B-35W	7.0x4.5	35
DSS A/B-35R	5.5x3.5	35
DSS A/B-45W	7.0x4.5	45
DSS A/B-45R	5.5x3.5	45
DSS A/B-55W	7.0x4.5	55
DSS A/B-55R	5.5x3.5	55



Disposable Umbilical Clamp and Cutter

DUC



Disposable Purse String Instrument

PSI

Code	Staple	Cartridge Length (mm)	Suture Length (mm)
PSI-45	18	47	600
PSI-65	28	67	700



Nasogastric Feeding Tubes



Code	Catheter			Guide Wire	
	Outer Circumference (mm)	Outer Diameter (mm)	Length (mm)	Diameter (mm)	Length (mm)
FTCH 5 A/B/C	5	1.67±0.5	510±80	0.5±0.2	≥60
FTCH 6 A/B/C	6	2.0±0.5	510±80		
FTCH 8 A/B/C	8	2.67±0.5	910±200	0.5±0.2	≥950
FTCH 10 A/B/C	10	3.33±0.5	910±200		
FTCH 12 A/B/C	12	4.0±0.8	1100±300		
FTCH 14 A/B/C	14	4.67±0.8	1400±300		

Plug meshes


Code	D(mm)	Height (mm)
SB-1	34	25
SB-2	36	26
SB-3	37.5	27
SB-4	39	29
SB-5	41.5	36
SB-6	40	35
SB-7	43	37
SB-8	45	38
SB-9	46.5	39
SB-10	48	41
SB-11	50	41

Flat meshes

Code	W (mm)	Height (mm)
25 x 100	25	100
50 x 100	50	100
50 x 110	50	110
50 x 120	50	120
60 x 110	60	110
60 x 150	60	150
70 x 110	70	110
80 x 110	80	110
80 x 150	80	150
100 x 100	100	100
100 x 150	100	150
150 x 150	150	150
200 x 200	200	200
200 x 250	200	250
200 x 300	200	300
250 x 300	250	300
300 x 300	300	300

Surgical Mesh and Hernia Plugs

Polypropylene Mesh (PP) Anti-adhesion mesh (PVDF+PP)





Disposable Hemorrhoid Ligator

DHL

DHL-A/B-15



Pelvic Organ Prolapse Repair System



Purse String with Needles

SZX

SZX-56/66 (Needle length 56mm, 66mm)


Disposable Circumcision Devices

Code	Inner Diameter (mm)	Tolerance (mm)
SHQQ-A/B 11	11	±1
SHQQ-A/B 12	12	±1
SHQQ-A/B 13	13	±1
SHQQ-A/B 14	14	±1
SHQQ-A/B 15	15	±2
SHQQ-A/B 16	16	±2
SHQQ-A/B 19	19	±2
SHQQ-A/B 20	20	±2
SHQQ-A/B 22	22	±2
SHQQ-A/B 24	24	±2
SHQQ-A/B 26	26	±2
SHQQ-A/B 28	28	±3



Disposable Circumcision Stapler

DCCS



Code	I.D (mm)	O.D (mm)	Staple	Height (mm)
DCCS-12	12	21	9	3.0
DCCS-15	15	24	11	3.0
DCCS-18	18	27	12	3.0
DCCS-20	20	29	14	3.0
DCCS-22	22	31	15	3.0
DCCS-24	24	33	16	3.0
DCCS-26	26	35	17	3.0
DCCS-28	28	37	19	3.0
DCCS-30	30	39	20	3.0
DCCS-32	32	41	22	3.0
DCCS-34	34	43	24	3.0
DCCS-36	36	45	25	3.0

Disposable Clip Applier

DCA (Ti Clip)

Code	Cannula Diameter (mm)	Metal Clip Specification (mm)	Close length of Metal Clip (mm)
DCA (A,B) 5-S	φ5	3.5×8.0	8.0
DCA (A,B) 10-L	φ10	5.0×9.0	9.0
DCA (A,B) 10-ML	φ10	4.5×9.0	9.0
DCA (A,B) 10-XL	φ10	4.5×11.0	11.0
DCA (A,B) 10-XXL	φ10	6.0×11.0	11.0

