

When ergonomic design  
and intuition come together



**ERAGON** *axial*

Axial handle unites with modular system

# ERAGON<sub>axial</sub>

When ergonomic design and intuition come together

The **ERAGON modular** forceps and scissor system, already having a proven track record, has been expanded by the addition of a completely new axial handle. **ERAGON axial** combines all the benefits of the existing handle concept, such as wide finger contact surfaces, optimum operation of the rotation knob, precise grip, and secure hold ratchet lock. **ERAGON axial** offers fatigue-free working and precision with a new level of handling. See for yourself.



## Ergonomic design gains upper hand

When the new handle was being created, simple handling and user-friendly operation were a top priority. The slightly angled contour permits a relaxed working position. The wide contact surfaces for the fingers facilitate fatigue-free working without the development of pressure points and therefore offer maximum surgeon comfort.



### Turning with finesse

The optional knob attachment made of silicone makes operation even easier. A convenient accessory that is helpful for all hand sizes – small and large alike. The reusable attachment is simple to remove and easily reprocessed.

### Safety under control

The intuitive, ratchet lock provides sensitive handling and a secure grip.

### Skillfully combined

- **ERAGON axial** is fully compatible with all the components in the **ERAGON modular Series**. All sheaths and jaw inserts can be combined as required.
- The **"Click-it"** system provides the intuitive assembly and disassembly familiar for ERAGON.
- Make the most of the multi-faceted diversity of the successful **ERAGON modular Series**.

# ERAGON<sub>axial</sub>

When ergonomic design and intuition come together



**ERAGON axial handle**  
with  
ratchet lock ..... 83930082



**ERAGON axial handle**  
without lock, with  
HF connector ..... 83930083



**ERAGON axial handle with**  
ratchet lock, with  
HF connector ..... 83930084



**ERAGON axial handle**  
with ratchet lock,  
not rotatable ..... 83930085



Option:  
**Knob attachment**  
**ERAGON axial,**  
autoclavable ..... 8988

## ERAGON<sub>modular</sub>

All ERAGON modular components  
are compatible with ERAGON axial.



WL = 450 mm

WL = 330 mm

WL = 260 mm

